

Patient: (42 years)
Recorded: Wednesday,
Mar 23, 2022, 8:35:09 AM
Heart Rate: 69 bpm
Duration: 30s

Kardia Advanced Determination:
Normal Sinus Rhythm

Kardia

Enhanced Filter, Mains Filter: 60 Hz Scale: 25mm/s, 10mm/mV




Patient: (42 years)
Recorded: Wednesday, Mar 23, 2022, 8:35:09 AM
Heart Rate: 69 bpm
Duration: 30s

Kardia Advanced Determination:
Normal Sinus Rhythm



R-R Interval Plot

Detection of QRS locations allows Kardia AI to measure the distance between heartbeats or the RR interval. This can be used to review variability of heart rate, which may be useful in understanding heart rate variability, or to visually display irregularity in rate (such as in Afib).



Patient: (42 years)
Recorded: Wednesday,
Mar 23, 2022, 8:35:09 AM
Heart Rate: 69 bpm
Duration: 30s

Kardia Advanced Determination:
Normal Sinus Rhythm



Average Beat Plot

The average beat display is the average of all the normal, non-ectopic (extra/missed) beats in the EKG. An average beat display is intended to be a simple visual representation of one beat in a 30-second EKG.

